


Introducing Beaty

For Effective CPR



Beaty

The power to save lives

Approximately 30,000 Australians die from Sudden Cardiac Arrest each year. Time is a crucial factor as early and effective CPR can dramatically increase the chance of survival and outcomes among victims.

Introducing Beaty

Beaty is an inexpensive, compact device that provides first responders with CPR feedback.

Smart Beaty provides real time, audible CPR feedback through measuring the depth of chest compressions during CPR.

Simple No buttons or displays and is small enough to fit in your pocket.

Affordable and priceless, you just might save a life!



How to use Beaty?



In case of CPR
call for HELP



Simply place Beaty in the
middle of the patient's chest



Push down until a
"beep" sound is heard

Beaty is a great addition to any first aid kit,
AED carry case or even on your key ring.

